

TREC 2022 Clinical Trials Track Runs

Run	Org	Run Type
ANIR_demo	CSIROmed	automatic
CSIROmedANIR	CSIROmed	automatic
doct5query	CSIROmed	automatic
monobert500	CSIROmed	automatic
zs_bert_500	CSIROmed	automatic
DoSSIER.1	DOSSIER	automatic
DoSSIER.2	DOSSIER	automatic
DoSSIER.3	DOSSIER	automatic
DoSSIER.4	DOSSIER	automatic
DoSSIER.5	DOSSIER	automatic
bm25_bi.filtered	els_dshs	automatic
bm25_st.bienc	els_dshs	automatic
senttr	els_dshs	automatic
st_distilbert	els_dshs	automatic
run1	ELUCIDATA_ML	automatic
frm3	h2oloo	automatic
frocchio	h2oloo	automatic
frocchio_monot5_e	h2oloo	automatic
frocchio_monot5_ed	h2oloo	automatic
rocchio	h2oloo	automatic
ims_BM25Filtered_kw	iiia-unipd	manual
ims_BM25Filtered_s	iiia-unipd	manual
ims_RM3Filtered_kw	iiia-unipd	manual
ims_RM3Filtered_s	iiia-unipd	manual
ims_T5summarizer	iiia-unipd	manual
jbnu1	jbnu	automatic
jbnu2	jbnu	automatic
phir1m1	phi_lab	automatic
phir2m2	phi_lab	automatic
phir3m1prf	phi_lab	automatic
phir4m1prf2	phi_lab	automatic
phir5m2prf	phi_lab	automatic
all_doc.search	POZNAN	automatic
condition.search	POZNAN	automatic
condition_title.search	POZNAN	automatic
IKR3_BSL	UNIMIB	automatic
IKR3_BSL_TT_HW	UNIMIB	automatic
IKR3_BSL_TT_MW	UNIMIB	automatic
IKR3_TT_BW	UNIMIB	automatic
IKR3_TT_MW	UNIMIB	automatic
crossrunner	V-Ryerson	manual